

# TPC

Ratchet Assembly teeth replacement procedure.

- Remove the two Allen screws shown.



- Once the Allen screws are removed, use a flat head screw driver and gently pry the top portion of the ratchet assembly apart.



# TPC

Ratchet Assembly teeth replacement procedure.

- Once the ratchet is apart, locate the two locking teeth.




- Gently lift up on each tooth and remove the spring fastener.



# TPC

Ratchet Assembly teeth replacement procedure.

- Once removed replace the locking teeth.



- Once both teeth are replaced, insert the locking lever . Rock the lever from left to right to secure it in place.



# TPC

Ratchet Assembly teeth replacement procedure.

- Insert ratchet stem back into assembly.



- Work the lever left and right to allow the gears to seat back into the assembly.



- Replace the two Allen screws that were removed in step one.