

# ***TPC***

## **Software Compatibility List**

The following imaging applications have been tested and will work with all of the AIC 5XXX intraoral cameras.

- Apteryx X-ray vision and X-ray Vision lite.
- Dexis Versions
- Patterson Eagle soft
- TPC Soft Image
- Professor Suni “Apteryx imaging”
- CliniView Quality Systems INC
- Dentimax “Apteryx Imaging”
- Kodak Imaging.
- Eva
- EZ 2000
- XDR Digital imaging.

&©j`¥

We have confirmed that the following applications WILL NOT integrate with the Advance cam Intraoral camera systems. “Freeze button will not perform a capture into the software application”.

- Visix Must purchase additional plug-in from software manufacture.
- ~~XXX~~ ¥³ ¥
- Tiger View

If you or your customer has an application that is not listed above please feel free to contact us to discuss the integration options.

This list will be updated as new information is brought to our attention. Some software updates may cause conflicts on how your camera operates. Please let us know if you experience any adverse issues with your camera if an update is causing issues with your integration.

### System requirements

- Windows XP, Vista and windows 7.
- 32 & 64 bit compatible hardware / software
- Minimum 512 MB ram. Only common on windows XP 32 bit OS.
- USB 2.0

There is currently no support for the advance cam on MAC / Linux OS.

TPC does not have any affiliation with any of the software vendors listed. This list is to be used as a reference only. If you have any questions please contact [service@tpcdental.com](mailto:service@tpcdental.com)