






Advance 500

TUBO CURING / BLEACHING LIGHT



- **R: RAMP CURING**
Ramp Curing allows versatile curing time setting in cycle of 20/40/60 seconds.
- **B: HIGH-SPEED CURING (BOOST FUNCTION)**
The light output with the boost function in using the 8-mm turbo tip will be 1,500 mw/cm² (average output); the curing time can be selected as 5 seconds or 10 seconds.
- **C: CONVENTIONAL CURING**
Conventional Curing function in cycle of 10/20/30/40/50/60 seconds.
- **BL: BLEACHING FUNCTION**
To be applied as bleaching function in cycle of 30 seconds curing mode. This function should be applied with suitable bleaching material.
- **BUILT-IN RADIOMETER**
- **LED DIGITAL DISPLAY**
- **LIGHT BULBS 12V, 75W**

LIGHTING GUIDES

	Straight 13mm diameter light guide for large facial surface
	Curved 13mm diameter 80° light guide for posteriorocclusal curing
	Curing 8mm diameter 60° light guide for general Curing
	Curing 3mm diameter 90° light guide for layer curing of posterior restoration
	Turbo light guide 8mm(tip) increase output up to 35% useful for deep or opaque restoration