


Model #800N



### Features

- Extremely gentle on patient, causes little or no sensitivity
- No heat build up in handpiece tip with extended use
- Autoclavable handpiece
- Extensive power range to allow widest choice of clinical application
- Fine water spray eliminates excessive overspray onto patient
- Thin subgingival tips available
- 3 tips included with unit
- 1 Year Warranty

### Additional Benefits

- The Advance 800N eliminates the need for external water line and power cord connections
- The Advance 800N operates from unit foot control - no additional foot control needed
- The Advance 800N allows for a clean appearing operatory eliminating all visible tubs power lines.

Shipping Dimensions: 11"L x 8"W x 7"H / 4 lbs.  
(16/cs.) 23"L x 15"W x 21"H / 58 lbs.

### Specifications

- Transformer: 24V, 1 amp
- Activated by unit foot control – no additional foot control needed
- Small electronic box, 2-3/8"L x 1.5"W x 1.25"H
- Easy tip system allows tips to be changed easily
- 5 tips included
- Absorbed power: max 15 watt
- Scaler: max 10 watt, 24-28 kHz, automatic tuning
- Water pressure: 30-145 PSI
- Water consumption: 10-50 cc/min
- No external water line / power cord

### Package includes the following:

- Advance 800N piezo ultrasonic scaler with handpiece
- #1, #2, #3, #4, #5 tips
- Tip wrench
- Transformer
- Control box

