

System Configurations



PC Setup Includes TPC Soft Image Software

AdvanceCAM USB Direct #AIC899

PC

PC Printer

Retail: \$1,390 (#AIC899)

TV / LCD Monitor Setup

AdvanceCAM #AIC888

"AV output"

#AIC750 Converter

TV or LCD

Retail: \$1,840 (#AIC888 and #AIC750)

PC Setup with Docking Station Includes TPC Soft Image Software

AdvanceCAM #AIC888

#AIC805 or #825 w/USB Port

PC

PC Printer

Retail: \$2,060 (#AIC888 and #AIC805)

TV / LCD Monitor Setup with Docking Station

AdvanceCAM #AIC888

AV or S-Video Output
VGA Output (#AIC825)

#AIC805 or #825 Docking Station

TV or LCD

Retail: \$2,060 (#AIC888 and #AIC805)

Wireless Network with PC

AdvanceCAM #AIC888 with #AIC900 Wireless Transmitter

#AIC815 or #AIC835 Wireless Docking Station with USB

PC

PC Server

Network PC Printer

PC

Retail: \$2,440 (#AIC888, #AIC900 and #AIC815)

Wireless Network with TV / LCD Monitor

AdvanceCAM #AIC888 with #AIC900 Wireless Transmitter

#AIC815 or #AIC835 Wireless Docking Station

AV or S-Video Output
VGA Output (#AIC835)

TV or LCD

Retail: \$2,440 (#AIC888, #AIC900 and #AIC815)

SD Memory Card Setup

AdvanceCAM #AIC888

#AIC650

AV or S-Video Output

TV or LCD

Photo Printer

*Place SD Card into Video Printer

Retail: \$2,080 (#AIC888, #AIC650)

Wireless SD Memory Card Setup

AdvanceCAM #AIC888 w/#AIC900 Wireless Transmitter

#AIC665

AV or S-Video Output

TV or LCD

Photo Printer

*Place SD Card into Video Printer

Retail: \$2,440 (#AIC888, #AIC900 and #AIC665)

SD Memory Card Device Setup *Place SD Card into Video Printer

#AIC600

Photo Printer

AV or S-Video Input

AV or S-Video Output

TV or LCD

Retail: \$720 (#AIC600)

Built-In Setup

AdvanceCAM #AIC888

AIC#300 Built-In Docking Station

#WIC840 Flat Screen Monitor Bracket

#MTFT15 Flat Screen Monitor

Retail: \$1,980 (#AIC888 and #AIC300)