

# LED 60


## Cordless Curing Light System

Model LED60



### Features

- Low battery signal – visual and audible alarm
- Easy to read – LCD display
- High light intensity – 5 to 10 second cure
- No down time
- Flexibility – corded and cordless light cure
- Quiet fan – extends life of LED
- Sleep mode – conserves power, saves battery life



### Specifications and Accessories

#### Switching Power Adapter

Power supply: AC 100-240V / 50-60 Hz  
 Dimensions: 10"L x 2"W x 8"H  
 Weight: 24 ounces

#### Hand Piece

Light source: 5W light emitting diode  
 Wavelength: 430-490 nm  
 Light Intensity: 1,400 nW/cm<sup>2</sup>  
 Battery: Lithium ion  
 Warranty: 1 year

### Standard Accessories

- LED 60 Handpiece
- Charging base
- 8mm light guide
- Eye protector
- Lithium ion battery
- Power supply
- User's Manual

Description	Model	Suggested Retail Price
LED 60 Cordless Curing Light	ALED-60	\$650

#### Package includes the following:

LED 60 Cordless Curing Light	✓
Charging Base	✓
Protective Eye-shield	✓
Countertop Stand with Power Indicator	✓
8mm Light Guide	✓
Built-in power transformer	✓
AC power supply	✓

### Options:

Options:	Retail
P108 Replacement light guide (8mm)	\$90.00
P208 Replacement protective eye-shield	\$16.00
ALED 60-BATT Replacement rechargeable battery for LED 60	\$100.00

Shipping Dimensions: 10"L x 7"W x 6"H / 3 lbs.  
 (8/cs.) 21"L x 14"W x 10"H / 16 lbs.