

# D650

## Digital Amalgamator



(D-650)

### Features

- ISO 9002 / ISO 13488 quality management certified
- Input voltage: 110 or 230V
- Touch pad selector
- Pre-programmed time settings
- Quiet, high-speed mix capsules
- Built-in automatic timer
- Accepts all capsules
- Customized trituration times
- Vibration frequency  $\geq$  4,400 RPMs
- 1 Year Warranty

### Description

D650 Digital Amalgamator

### Model

D-650

### Suggested Retail Price

\$390

### Package includes the following:

D650 digital amalgamator  
AC power cord



Shipping Dimensions: 12”L x 12”W x 12”H / 8 lbs.