



Model #L700

### **Features:**

- Light may be adjusted diagonally, vertically, and horizontally
- Sharply defined rectangular light pattern
- 2-Position intensity switch
- Dichroic coated reflector filters out infrared rays providing cool light
- Easy open safety shield allows aseptic cleaning
- Removable splash shield
- Color correct
- 1 year warranty
- Packing dimensions: 26" x 14" x 8" / 20 lbs.

### **SPECIFICATIONS**

Input Voltage:	110 V, 60Hz
Output Voltage:	12 V High, 10.5 V Low
Power Output:	55 W High, 50 W Low
Bulb Type:	Tungsten Halogen Type 12 V - 55 W
Light Output:	24,000 Lux High 12,000 Lux Low
Color Temperature:	4500° Kelvin
Focal Distance:	27"
Light Pattern:	3" x 7" @ 27"

