






Model LED 39N

Features

- Lightweight
- Three modes of operation: fast, ramp, pulse
- Ergonomic contoured design



Specifications and Accessories

Switching Power Adapter

Power supply: AC 100-240V / 50-60 Hz
 Dimensions: 8"L x 1"W x 1"H
 Weight: 5 ounces
 Warranty: 2 years

Handpiece

Light source: 5W light emitting diode
 Wavelength: 420-480 nm
 Light Intensity: 1,100 mW/cm²
 Warranty: 2 years



Description

Description	Model
LED 39N Curing Light	LED 39N

Package includes the following:

LED 39N Curing Light	✓
Stand	✓
8mm Light Guide	✓
Eye Protector	✓
Power Supply	✓
User's Manual	✓

Options:

P109	Replacement light guide (8mm)
P219	Replacement protective eye shield
LED39N-PS	Replacement power supply
LM200	Light Meter

Standard Accessories

- LED 39N Handpiece
- Stand
- 8mm light guide
- Eye protector
- Power supply
- User's Manual

Shipping Dimensions: 9"L x 4"W x 2.5"H
 2 lbs.

