



Model LED 60N



### Specifications and Accessories

#### Switching Power Adapter

Power supply: AC 100-240V / 50-60 Hz  
 Dimensions: 10"L x 2"W x 8"H  
 Weight: 24 ounces

#### Handpiece

Light source: 5W light emitting diode  
 Wavelength: 430-490 nm  
 Light Intensity: 1,400 mW/cm<sup>2</sup>  
 Battery: Lithium ion  
 Warranty: 2 years



### Standard Accessories

- LED 60N Handpiece
- Charging base
- 8mm light guide
- Eye protector
- Lithium ion battery
- Power supply
- User's Manual

Shipping Dimensions: 10"L x 7"W x 6"H / 3 lbs.  
 (8/cs.) 21"L x 14"W x 10"H / 16 lbs.

