

Assistant's "North Carolina" Mobile Cabinet

#TMC-318-120



- 21.5" W x 19" D x 32" H
- solid top with pull handle
- amalgamator well drop leaf door
- 3 instrument drawers / 1 storage drawer
- 1 year warranty
- color: lite grey or white
- packing dimensions: 40" x 40" x 36", 110 lbs.

Assistant's "Alabama" Mobile Cabinet

#TMC-317-160



- 21.5" W x 19" D x 32" H
- sliding counter top with amalgamator well
- 1 instrument drawer / 2 storage drawers
- side to side sliding top
- 1 year warranty
- color: lite grey or white
- packing dimensions: 40" x 40" x 36", 110 lbs.

Utility Mobile Cabinet

#TMC-653-180



- 14.5" W x 16" D x 32" H
- Designed for ultrasonic scalers, curing lights, electrosurge, intraoral cameras, etc.
- 1 instrument drawer / 2 adjustable shelves
- 1 year warranty
- Color: lite grey or white
- Packing dimensions: 40" x 40" x 36", 65 lbs.

Doctor's Mobile Cabinet

#TMC-310-140



- 21.5" W x 19" D x 32" H
- 4 instrument drawers / 3 storage drawers
- 1 year warranty
- Color: lite grey or white
- Packing dimensions: 40" x 40" x 36", 110 lbs.

