



15" Deluxe Flat Screen Monitor

WITH BUILT IN X-RAY VIEWER

Active-Matrix TFT LCD Monitor


- 1074 x 768 resolution
- 500:1 high contrast ratio
- 160-degree viewing angle

Mirror Image Function

- Enables viewing from opposite perspective
- Better communication between Doctor and patient

Built-In Stereo Sound System

- Dual built-in stereo speakers



(#MTFT15)

Anti-Flash and Anti-Glisten Intensified Glass Sheet

- Protects thin film transistor LCD panel
- Longer LCD life

Signal Hub Device

- Monitor for: pc, satellite tv, cable tv, dvd, vcd, digital camera, video camera, intraoral camera

Wall, Post or Counter Top Mounted

- Multiple monitors may be networked
- VESA wall mountable

Built-In Buffer Storage for Images

- Built-in storage for immediate comparison before and after treatment



Built-In Xray Viewer

- Dry or wet film
- Enhanced 50:1 enlargement
- Images easily stored into computer



Versatile One-Touch Remote Control

- Built-in infrared remote mouse, decimal channel selector, laser pointer, mirror imaging, function selector, adjustment for contrast, luminance and color.



(#MTFT15 with optional #WIC840)