

Model PM25



Features

- 25Khz
- Precision automatic tuning circuitry
- Linear power control for exact power setting
- Accepts Cavitron®, P, SLI, TFI, and FSI scaler inserts
- Compact and lightweight
- Auto-ranging 110-240 VAC input
- Includes one ultrasonic scaler insert and foot control
- 1 year warranty

Power you want,
where you want it!

Specifications

Weight: 6 lbs.
Dimensions: 8"(W) x 8"(D) x 3"(H)
Output power to tip: Operatory adjust approximately 10-40 watts
Handpiece cable: 7.5'
Power supply to unit: 110-240 VAC Input, 50/60 Hz
Water flow and supply pressure: 18 cc to 30 cc per minute flow
from insert tip
Supply pressure to scaler is 25-60 PSIG city water or bottle system

