

#DR-1102  
Doctor's Stool



### Features

- Three seat function
- Seat pan tilt
- Six-way lumbar support
- Two height ranges
- Height, back and seat adjustment levers
- Durable tile/carpet casters

### Dimensions

Backrest vertical adjustment range: 0" to 3.25"

Seat width: 18"

Seat height range: 18.25" to 25"

5-leg base diameter: 24"

Packing dimensions: 24" x 22" x 15"

Weight: 35 lbs.

#DR-5109  
Doctor's Stool with Armrest



### Features

- Three seat function
- Adjustable armrest
- Seat pan tilt
- Six-way lumbar support
- Two height ranges
- Height, back & seat adjustment levers
- Durable tile/carpet casters

### Dimensions

Backrest vertical adjustment range: 0" to 3.25"

Seat width: 18"

Seat height range: 18.25" to 25"

5-leg base diameter: 24"

Packing dimensions: 24" x 22" x 15"

Weight: 35 lbs.

#AS-1101  
Assistant's Stool



### Features

- Six-way lumbar support
- Height, back and seat adjustment levers
- Adjustable foot ring
- Durable tile/carpet casters

### Dimensions

Abdominal support vertical adjustment range: 0" to 3.25"

Seat width: 18"

Seat height range: 21.5" to 27.5"

5-leg base diameter: 24"

Packing dimensions: 24" x 22" x 15"

Weight: 45 lbs.

