



UL110

3 Curing Modes are Available:

- 1- Fast (F): Full light output
- 2- Ramp (R): Soft start light output
- 3- Pulse (P): On/Off light every 0.2 seconds

- Rapid Curing with high light intensity
- No heat curing: LED never heats up
- No noisy fan: wonderful quiet curing
- Long use-life LED: 100,000 hours
- Uni LED may be built-in to dental unit
- Spectrum range: 450-490nm
- Light weight: 76g (without light guide)
- Switch power supply: free use from 100V to 240V
- 1 year warranty

SPECIFICATIONS & ACCESSORIES

Switching Power Adapter

Input: AC 100-240, 47,63 Hz
Output: DC 5V, 1.5A
Output Power: 7.5W Max

Light Source

Model Type: High-flux Blue LED
Wavelength: 430-490 nm

Curing Handpiece

Length: 10.5"
Diameter: 1"
Net weight: 94.5g

