



#S1025 – 25K Universal Scaler Insert  
(metal handle)

#S1030 – 30K Universal Scaler Insert  
(metal handle)

### Features

- External water
- Interchangeable for use with other brands of scalers
- 90 day warranty against defects arising from faulty materials or workmanship



#S2025 – 25K Universal Slim Scaler Insert  
(metal handle)

#S2030 – 30K Universal Slim Scaler Insert  
(metal handle)

### Features

- Slim tip design
- External water
- Interchangeable for use with other brands of scalers
- 90 day warranty against defects arising from faulty materials or workmanship



#S3025 – 25K Universal Scaler Insert  
(plastic handle)

#S3030 – 30K Universal Scaler Insert  
(plastic handle)

### Features

- Internal water delivery system delivers water directly to the insert tip and focuses it where it is needed most
- Large, round, non-slip plastic grip for greater tactile sensitivity, rotational control and reduced hand fatigue
- Better water management and less water usage means greater patient comfort
- Interchangeable for use with other brands of scalers
- 90 day warranty against defects arising from faulty materials or workmanship



#S4025 – 25K Universal Slim Scaler Insert  
(plastic handle)

#S4030 – 30K Universal Slim Scaler Insert  
(plastic handle)

### Features

- Internal water delivery system delivers water directly to the insert tip and focuses it where it is needed most
- Large, round, non-slip plastic grip for greater tactile sensitivity, rotational control and reduced hand fatigue
- Better water management and less water usage means greater patient comfort
- Interchangeable for use with other brands of scalers
- 90 day warranty against defects arising from faulty materials or workmanship

