




### Features

- Light head will rotate to any position
- Quiet fan operation
- Multi-arch whitening application
- Four high intensity LED's
- Digital display with audio feedback
- High power blue light emitting diodes
- Compact mobile stand for easy storage
- Adjustable bleaching time with microprocessor-controlled timer
- Smart auto-select power: 110-240 Volts AC, 50/60 Hz
- LED bleaching system includes illuminator, mobile base, adjustable arm, power cord & user manual
- 1 year warranty
- Packing dimensions: 22" x 9" x 10"
- Weight: 15 lbs.

### Specifications

Broad Spectrum .....	430-490 nm & 395-415 nm
Light Intensity .....	1,500 mW/cm <sup>2</sup>
Voltage Output .....	110-240 V, 50/60 Hz
Dimensions .....	60" (H) x 20" (W)
Weight .....	15 lbs.



Model #2505