



**Winter 2025**

**INNOVATION•RELIABILITY•QUALITY**

SOLD THROUGH AUTHORIZED DEALERS ONLY

Prices in **RED** denote Show Special pricing

Show Specials **25% OFF**



**50% OFF**  
Factory Installed  
Built in Electric Motor  
Retail: \$1,790  
**\$895 NOW!!**

## Bellagio Operatory Package

#BP5000-L7601LED-2.0

Retail: \$17,180

Show Special **\$12,885**

#BP5015-L7601LED-2.0

Retail: \$16,480

Show Special **\$12,360**



**50% OFF**  
Factory Installed  
Built in Electric Motor  
Retail: \$1,790  
**\$895 NOW!!**

## Mirage Operatory Package

#MP2015-600LED-2.0

Retail: \$13,080

Show Special **\$9,810**

#MP2000-600LED-2.0

Retail: \$13,780

Show Special **\$10,335**



**50% OFF**  
Factory Installed  
Built in Electric Motor  
Retail: \$1,790  
**\$895 NOW!!**

## Laguna Operatory Package

#LP2115-L550LED-2.0

Retail: \$10,080

Show Special **\$7,560**

#LP2100-L550LED-2.0

Retail: \$10,780

Show Special **\$8,085**



**50% OFF**  
Factory Installed  
Built in Electric Motor  
Retail: \$1,790  
**\$895 NOW!!**

## Mirage Swing Operatory

#MSP3500-L600LED-2.0

Retail: 14,456

Show Special **\$10,842**

**TPC**  
5/2 Year  
WARRANTY



**50% OFF**  
Factory Installed  
Built In Electric Motor  
Retail: \$1,790  
**\$895 NOW!!**

### *Bellagio Operatory Package*

#BP5000-L700LED-2.0

**Retail: \$15,570**

**Show Special \$11,680**

#BP5015-L700LED-2.0

**Retail: \$16,270**

**Show Special \$12,205**

**TPC**  
5/2 Year  
WARRANTY



**50% OFF**  
Factory Installed  
Built In Electric Motor  
Retail: \$1,790  
**\$895 NOW!!**

### *Mirage Operatory Package*

#MP2015-L700LED-2.0

**Retail: \$13,570**

**Show Special \$10,180**

#MP2000-L700LED-2.0

**Retail: \$14,270**

**Show Special \$10,705**

**TPC**  
5/2 Year  
WARRANTY



Optional:  
LCD monitor  
Monitor bracket

**50% OFF**  
Factory Installed  
Built In Electric Motor  
Retail: \$1,790  
**\$895 NOW!!**

### *Laguna Operatory Package*

#LP2115-L700LED-2.0

**Retail: \$10,770**

**Show Special \$8,080**

#LP2100-L700LED-2.0

**Retail: \$11,470**

**Show Special \$8,605**

**TPC**  
5/2 Year  
WARRANTY



**50% OFF**  
Factory Installed  
Built In Electric Motor  
Retail: \$1,790  
**\$895 NOW!!**

### *Mirage Swing Operatory*

#MSP3500-2.0

**Retail: 12,570**

**Show Special \$9,430**



**#4000-2.0-  
Mirage 2.0 Hydraulic Chair  
Cast Base**

Retail: \$8,590  
**Special \$6,442**

**#3000-2.0-  
Mirage 2.0 Hydraulic Chair  
Cast Base (No Armrest)**

Retail: \$8,390  
**Special \$6,293**



**#4000-2.0  
Mirage 2.0 Hydraulic Chair**

Retail: \$7,390  
**Special \$5,543**

**#3000-2.0  
Mirage 2.0 Hydraulic Chair  
(No Armrest)**

Retail: \$7,190  
**Special \$5,393**



**#L2000-2.0  
Laguna 2.0 Electromechanical  
Patient Chair**

Retail: \$4,590  
**Special \$3,442**

**#L3000-2.0  
Laguna 2.0 Electromechanical  
Patient Chair (no armrests)**

Retail: \$4,490  
**Special \$3,367**



**#LOP2000-L550LED-2.0  
Laguna 2.0 Orthodontic Package**

Retail: \$9,416

**Special \$7,062**



**#PC2630  
Portable Dental Unit**  
• 2 hp manual selection  
• 3 way syringe  
• Water bottle system  
• Waste container  
• 3/4 hp compressor  
• Air venturi vacuum

Retail: \$2,390  
**Show: \$1,792**



**#PC2930  
Portable Dental Unit**  
• 2 hp automatic selection  
• **LED 39 curing light**  
• **Built in piezo scaler**  
• 1 year warranty

Retail: \$3,300  
**Show: \$2,475**



**#PC2750  
Tripod LED Light**  
• LED illumination  
• Tripod with carry bag  
• Intensity adjustment  
• 2 year warranty

Retail: \$378  
**Show: \$284**



**#PC2700  
Portable Dental Chair Package**

Model	Retail	Show
• PC 2720 Patient Chair	\$1,090	<b>\$820</b>
• PC 2740 Dr. Stool	\$480	<b>\$360</b>
• PC 2760 LED Light	\$260	<b>\$195</b>
• PC 2780 Tray / Spittoon	\$290	<b>\$218</b>
• PC 2790 Carry Bag	\$190	<b>\$143</b>
• PC 2795 Stool Carry Bad	\$100	<b>\$75</b>

### #L7601-LED

#### Luminous LED Light

- Post mounted
- LED illumination
- Lux: 3,000 - 35,000
- Intensity adjustment
- Toggle switch on / off
- 2 year warranty

Retail: \$2,900

**Show Special \$2,175**




### #L7605-LED

#### Luminous LED Light

- Sensor on / off

Retail: \$3,190

**Show Special \$2,395**



### #L7901-LED

#### Luminous LED Ceiling mount

Retail: \$4,390 **Show Special \$3,295**

### #R6101-LED

#### Radiant LED Operatory Light

- Post mounted
- LED illumination
- Lux: 3,000 - 35,000
- Intensity adjustment
- Toggle switch on / off
- 2 year warranty

Retail: \$1,990

**Show Special \$1,495**



### #R6105-LED

#### Radiant LED Operatory Light

- Sensor on / off

Retail: \$2,190

**Show Special \$1,645**

### #L600-LED

#### Lucent LED Operatory Light

- Post mounted
- LED illumination
- Lux: 3,000 - 35,000
- Intensity adjustment
- Toggle on / off
- 2 year warranty

Retail: \$1,500

**Show Special \$1,125**



### #L690-LED

#### Lucent LED Operatory Light

- Ceiling mounted

Retail: \$2,160

**Show Special \$1,620**



### #L700-LED

#### Vivid LED Operatory Light

- Post mounted
- LED illumination
- Lux: 3,000 - 35,000
- Intensity adjustment
- Sensor +Toggle on / off
- Composite cure mode
- 2 year warranty

Retail: \$1,590

**Show Special \$1,193**



### #L795-LED

#### Vivid LED Operatory Light

- Ceiling mounted

Retail: \$2,250

**Show Special \$1,680**



### #L625-LED

#### Lucent LED Operatory Light

- Post mounted
- LED illumination
- Lux: 3,000 - 35,000
- Intensity adjustment
- Sensor on / off
- 2 year warranty

Retail: \$1,700

**Show Special \$1,275**



### #L695-LED

#### Lucent LED Operatory Light

- Ceiling mounted

Retail: \$2,360

**Show Special \$1,770**



### #L550-LED

#### Lustrous LED Operatory Light

- Post mounted
- LED illumination
- Lux: 8,000 - 30,000
- Intensity adjustment
- Toggle / sensor
- 2 year warranty

Retail: \$1,300

**Show Special \$1,195**



### #L570-LED

#### Lustrous LED Operatory Light

- Ceiling mounted

Retail: \$1,960

**Show Special \$1,470**



### #5000-2.0

#### Bellagio Delivery System

- 3 hp position w Touch pad
- Cuspidor
- Assistants Package w Touch pad
- Hve / SE
- Wet dry foot control
- Built in Piezo scaler
- Built in LED curing light
- 3 way syringe
- Water bottle system
- 5 year warranty



**Retail: \$5,690**

**Show Special \$4,270**

### #5015-2.0

#### Bellagio Delivery System (no cuspidor)

**Retail: \$4,990**

**Show Special \$3,745**

### #MC500-2.0

#### Mirage Mobile Cart

- 3 hp position
- 3 way syringe
- Wet dry foot control
- Water bottle system
- Utility center with j-box
- 5 year warranty



**Retail: \$2,500**

**Show Special \$1,875**

### #MC501-2.0

#### Mirage Mobile Cart

- 3 hp position
- 3 way syringe, HVE + SE
- Wet dry foot control
- Water bottle system
- 5 year warranty

**Retail: \$2,800**


**Show Special \$2,100**

### #MS3000

#### Mirage Side Delivery

##### Cabinet Mount

- 3 hp position
- 3 way syringe
- Wet dry foot control
- Water bottle system
- Utility center
- 5 year warranty



**Retail: \$3,300**

**Show Special \$2,475**

### #2015-2.0

#### Mirage 2.0 Delivery System

- 3 hp position w touch pad
- Assistants Package w touch pad
- HVE / SE
- Wet dry foot control
- 3 way syringe
- Water bottle system
- 5 year warranty



**Retail: \$4,190**

**Show Special \$3,142**

### #2000-2.0

#### Mirage 2.0 Delivery System (with cuspidor)

**Retail: \$4,890**

**Show Special \$3,667**

### #MR4500-2.0

#### Mirage 2.0 Rear Delivery with Vacuum Package

- 3 hp position
- 3 way syringe
- Wet dry foot control
- Water bottle system
- Utility center
- 5 year warranty
- HVE
- SE
- Vacuum trap



**Retail: \$2,990**

**Show Special \$2,242**

### #MR4000-2.0

#### Mirage 2.0 Rear Delivery

- 3 hp position
- 3 way syringe
- Wet dry foot control
- Water bottle system
- Utility center
- 5 year warranty



**Retail: \$2,200**

**Show Special \$1,650**

### #OMC2375 Orthodontic Delivery

- Three way syringe
- HVE / SE valves
- Wet / dry foot control
- Water bottle system
- Utility center junction box
- 5 Storage drawers
- 1 year warranty
- Requires central air and vacuum



**Retail: \$3,190**

**Show Special \$2,390**



OMC-ST81



OMC-ST82



OMC-ST85



OMC-ST86

**Retail: \$420**

**Show Special \$315**

### #OMC2375CV Orthodontic Delivery

- Three way syringe
- HVE / SE valves
- Wet / dry foot control
- Water bottle system
- Waste bottle system w/ solids collector
- 5 Storage drawers
- 3/4 hp oil-less motor
- 1 year warranty
- Self contained



**Retail: \$4,700**

**Show Special \$3,525**

### #OMC2375CV-SL Orthodontic Delivery

- Built in Piezo scaler
- Built in LED 39 curing light



**Retail: \$5,600**

**Show Special \$4,200**

### #MC-3600CV

#### Mirage Mobile Delivery Cart

- Self contained
- Three way syringe
- Wet / dry foot control
- Water bottle system
- Waste bottle system
- Hve / SE valves
- 1 year warranty



**Retail: \$4,270**

**Show Special \$3,203**

### #MC-3600CV-SL

#### Mirage Mobile Delivery Cart

- Self contained
- Three way syringe
- Wet / dry foot control
- Water bottle system
- Waste bottle system
- Hve / SE valves
- Built in Piezo scaler
- Built in LED 39 curing light
- 1 year warranty



**Retail: \$5,170**

**Show Special \$3,878**

### #PC2635

#### Portable Delivery System

- Three way syringe
- Water bottle system
- 2 HP positions
- Disc style foot control
- Carry handle
- Pull handle
- 1 year warranty



**Retail: \$1,790**

**Show Special \$1,343**

### #PC2530

#### Mobile Vacuum System

- 2 HVE valves
- 1 SE valves
- solids collector trap 2 1/8"
- 7/8 hp vacuum motor
- 2L Waste container
- 4 large durable casters
- 1 year warranty



**Retail: \$2,100**

**Show Special \$1,575**



- #DR-7000  
Solace Doctor's Stool**
- 140 mm gas cylinder
  - Backrest tilt adjustment
  - Wide contoured seat pan
  - 5 dual wheel casters
  - 1 year warranty

**Retail: \$818**

**Show Special \$615**



- #DR-1102  
Mirage Doctor's Stool**
- 140mm gas cylinder
  - Backrest tilt adjustment
  - Wide seat pan
  - 5 dual wheel casters
  - 1 year warranty

**Retail: \$650**

**Show Special \$489**



- #HY-5107  
Doctor's Stool**
- 140 mm gas cylinder
  - Backrest tilt adjustment
  - Triangle seat pan
  - 5 dual wheel casters
  - 1 year warranty

**Retail: \$690**

**Show Special \$518**



- #SS-5105  
Mirage Saddle Stool**
- 200 mm gas cylinder
  - 5 dual wheel casters
  - 1 year warranty

**Retail: \$690**

**Show Special \$518**



- #DR-5045  
Laguna Doctor's Stool**
- 140mm gas cylinder
  - Backrest tilt adjustment
  - Wide seat pan
  - 5 dual wheel casters
  - 1 year warranty

**Retail: \$480**

**Show Special \$360**



- #AS-1101  
Mirage Assistant's Stool**
- 200 mm gas cylinder
  - Backrest tilt adjustment
  - Wide seat pan
  - Body support bar
  - 5 dual wheel casters
  - 1 year warranty

**Retail: \$818**

**Show Special \$615**



- #AS-5047  
Laguna Assistant's Stool**
- 200 mm gas cylinder
  - Body support bar
  - Circular seat pan
  - 5 dual wheel casters
  - 1 year warranty

**Retail: \$590**

**Show Special \$443**




- #XR-6101  
X-Ray Chair**
- Adjustable headrest
  - Pneumatic pump for base
  - Arm rest
  - Leg rest
  - 1 year warranty

**Retail: \$1,550**

**Show Special \$1,165**


**Standard Vinyl Upholstery Color Selection**



**Morbern Allante FFree Standard Vinyl**



**Optional NaguaSoft Vinyl Upholstery Upgrade:**



**Dr. Stool**

Retail \$148

**\$111**

**Asst. Stool**

Retail \$170

**\$128**

## Mechanical Room Combo Bundle



DC8228



WV3

SuperbAIR #DC8124 2HP Single Air Compressor (3 User)

OR

SuperbAIR #DC8228 4HP (Twin) Air Compressor (7 User)

+

Any model Superb VAC Vacuum Pump



DC8124



WV10

# 30% OFF Retail

## Mechanical Room

### TPC Superb Vac Wet Ring Vacuum System

Model		Show Special
WV3	(Single) 1 HP 110V or 220V	\$2,395
WV3R	(Single) 1 HP 220V w/ Recycler	\$2,695
WV5	(Single) 2 HP 220V	\$2,695
WV5R	(Single) 2 HP 220V w/ Recycler	\$2,995
WV6	(Twin) 1 HP 220V	\$4,945
WV6R	(Twin) 1 HP 220V w/ Recycler	\$5,245
WV10	(Single) 2 HP 220V	\$5,545
WV10R	(Single) 2 HP 220V w/ Recycler	\$5,845



### TPC Superb Air Oil-Less Air Compressor

Model		Show Special
DC8112	(Single), 1 HP 110V (1 Users)	\$4,465
DC8124	(Single), 2 HP 220V (3 Users)	\$4,945
DC8224	(Twin), 2 HP 220V (3 Users)	\$5,995
DC8226	(Twin), 3 HP 220V (5 Users)	\$7,195
DC8228	(Twin), 4 HP 220V (7 Users)	\$8,695



## Mobile Cabinet Systems



**#TMC-120 Assistant's North Carolina Mobile Cabinet**  
 • 21.5"W x 19"D x 32"H  
 • Color: White or light grey  
 • Solid top with pull handle  
 • Amalgamator well with drop leaf door  
 • 3 instrument drawers  
 • 1 storage drawer  
 • 1 year warranty  
 Shipping dimensions: 25" x 23" x 36" 112 lbs  
**Retail Price: \$990**  
**Show: \$745**



**#TMC-140 Doctor's Mobile Cabinet**  
 • 21.5"W x 19"D x 32"H  
 • Color: White or light Grey  
 • Solid top with pull handle  
 • 4 instrument drawers  
 • 3 storage drawer  
 • 1 year warranty  
 Shipping dimensions: 27" x 25" x 36" 116 lbs  
**Retail Price: \$990**  
**Show: \$745**



**#TMC-160 Assistant's Mobile Cabinet**  
 • 21.5"W x 19"D x 32"H  
 • Color: White or light Grey  
 • Solid top with pull handle  
 • Sliding counter top with amalgamator well  
 • 1 instrument drawers  
 • 2 storage drawer  
 • 1 year warranty  
 Shipping dimensions: 25" x 23" x 35" 115 lbs  
**Retail Price: \$990**  
**Show: \$745**



**#TMC-180 Utility Mobile Cabinet**  
 • 14.5"W x 16"D x 32"H  
 • Color: White or light Grey  
 • 1 instrument drawers  
 • 2 Adjustable shelves  
 • 1 year warranty  
 Shipping dimensions: 24" x 24" x 37" 72 lbs  
**Retail Price: \$478**  
**Show: \$360**