

# TPC

# Mirage 2.0

Chair Mounted Delivery System

**TPC  
5 Year  
Limited  
Warranty**

Model #2000



### Packing Dimensions

Unit:	51" x 21" x 17"	60 lbs.
PMU:	23" x 19" x 21"	36 lbs.
Cuspidor:	18" x 14" x 14"	20 lbs.



### Cuspidor Bowl:

The deep contoured bowl is made of fine porcelain for easy cleaning and infection control. The cuspidor operates smoothly with the timed bowl rinse, cup filler and swing-away arm.



### Optional:

Multi media flat screen monitor for chair side intraoral imaging, patient education and DVD viewing.



### Control Head:

Easy access controls for simple adjustments. Ultra soft handpiece tubing minimizing hand and wrist stress. 4 handpiece automatic air diaphragm block

**Standard Features**

- Asepsis automatic control for three handpieces
- Individual handpiece pressure adjustments
- Non-retracting water coolant with on/off
- Individual water coolant flow controls for all handpieces
- Air coolant flow control
- Handpiece tubing flush system
- Ultra soft handpiece tubing
- 6-pin straight gray asepsis tubing with Midwest connector
- Silicone non-skid pad
- Disc foot control wet/dry toggle with 7' of smooth tubing
- Equipment utility center with asepsis housing
- Vacuum canister with solids collector screen
- Water outlet quick disconnect
- Air outlet quick disconnect for prophy jet
- Self contained water system with city/bottle selector switch
- Straight light post for unit mount dental light
- Pre-wiring for fiber optic light source
- Air-actuated counterbalanced arm lock
- 3 Way Syringe on unit head
- Low voltage terminal block

# Mirage 2.0

**Chair Mounted Delivery System**



Model #2015



**Gravity Drain Cuspidor**

- Contoured bowl is made of fine porcelain for easy cleaning and infection control. The cuspidor operates with adjustable timed bowl rinse and cup filler
- Swing away cuspidor bowl

**Assistant's Instrumentation on Breakaway**

**Telescoping Arm**

- 1 autoclavable HVE with asepsis tubing
- 1 autoclavable saliva ejector tubing

**Options**

- Mounting bracket to mount flat screen monitor on light post
- Super thin X-ray viewer to mount on unit control head
- Built-in Led Curing Light
- Built-in Intraoral Camera
- Built-in Scaler
- Illumination Light Source
- Prophy Jet
- Post-mounted plastic utility tray
- Unit head mounted touch pad control

**Packing Dimensions**

Unit:	51" x 21" x 17"	60 lbs.
PMU:	23" x 19" x 21"	36 lbs.



Model #2000



**Standard:**  
Built in self contained 1500m water bottle system including city/bottle selector switch.



**Standard:**  
Built-in water quick disconnect for ultrasonic scaling