



DC6112D/A



DC6224D/A



# Superb Air Oil-less Air Compressor

The TPC Superb Air Oil-less Air Compressors meet the high standards of the dental professionals, providing the reliability you need from a brand you can trust. The new motor design delivers a reduced noise with a 5-year warranty. The **desiccant air drying system** efficiently filters out compressed air contaminants and moisture from the compressed air. Purge tank removes residual humidity in the compressor following each cycle.

Model	DC6112D/A	DC6122D/A	DC6214/DA	DC6224D/A	DC6326D/A	DC6228D/A
User Capacity	1-3users	1-3users	3-4users	5users	6-7users	8-9users
CDFM@80psi	4.9	4.9	9.8	9.8	14.7	19.6
# of Heads	1	1	2	2	3	2
Total HP	1.15	1.15	2.30	2.30	3.45	4.80
Voltage	110V	230V	110V	230V	230V	230V
Amps	<b>8</b>	<b>4</b>	<b>16</b>	<b>8</b>	<b>12</b>	<b>16</b>
Breaker Amps	20	20	30	20	30	30
dBA	62	62	62	62	72	72
Tank Capacity	10 Gallon	10 Gallon	16 Gallon	16 Gallon	23 Gallon	30 Gallon
Size (HxWxD)	16x23x27	16x23x27	26x21x29	26x21x29	34x22x29	41x22x36
Shipping Dimension (WxDxH)	29x23x35	29x23x35	34x26x35	34x26x35	39x25x38	46x27x41

### Features:

- New motor design to reduce noise and increase reliability
- Quiet operation
- Dry air prevents bacteria growth and air if filtered to 0.01 micron keeping 99.9997% of air contaminants out of your air system
- Energy efficient thermally protected motors
- Superb Air Oil-less compressors are available in single or twin heads
- **Desiccant air drying system to remove moisture from compressed air**
- **5-year warranty**

### Description

### Model

Superb Air Compressor, 1.15 HP, 110V	DC6112
Superb Air Compressor, 1.15 HP, 230V	DC6122
Superb Air Compressor, 2 x 1.15HP, (2.3HP) 110V	DC6214
Superb Air Compressor, 2 x 1.15HP, (2.3HP) 230V	DC6224
Superb Air Compressor, 3 x 1.15HP (3.45HP)230V	DC6326
Superb Air Compressor, 2 x 2.4 HP, (4.8HP)230V	DC6228

### Options

Control Panel – 2 Button- 1 Vac / 1 Air	#CP2AV
Control Panel – 3 Button- 1 Vac / 1Air / 1 Water	#CP3
Control Panel – 3 Button- 2 Vac / 1 Air	#CP3AV
Control Panel – 4 Button- 2 Vac / 1 Air / 1 Water	#CP4
Air Compressor and Vacuum Pump Stacking Rack	#CV22

### Maintenance kits for

	Model
DC 6112 / 6122 (air filter, desiccant filter, 5 micron filter)	DC6-PM1
DC 6214 / 6224 (air filter, desiccant filter, 5 micron filter)	DC6-PM2
DC 6226 (air filter, desiccant filter, 5 micron filter)	DC6-PM3
DC 6228 (air filter, desiccant filter, 5 micron filter)	DC6-PM4

\*SuperbAIR air compressors not to be used for CAD / 3D printers\*

### Shipping Dimensions: L X W X H

DC6112: 29" X 23" X 35"
DC6122: 29" X 23" X 35"
DC6214: 34" X 26" X 35"
DC6224: 34" X 26" X 35"
DC6326: 39" X 25" X 38"
DC6228: 46" X 27" X 41"