



Infinity 5 Intraoral Scanner Package Contents Guide





Revision Record

Version	Modify date	Modify content
1.1	2025.9.29	Added descriptions for the battery and power adapter



Product Name: Intraoral Scanner

Model: Infinity 5


Packaging unit: 1 set

1 Components

N o.	Component name	Specificati ons	Quantity	Pictures
1	Scanner Handpiece	Infinity 5	1	
2	Scanner Tip	/	4 (3standar d 1mini)	 
3	Protective Tip	/	1	
4	Base	/	1	
5	Battery	/	3	

6	Charging Holder		1	
7	Power Adapter		1	
8	Wi-Fi adapter		1	
9	Certificate Warranty Card	/	1	
10	U Disk with Software	/	1	
12	Packing box	/	1	
13	Practice model	/	1	

2 Label

No.	Label	Location
1		
2		
3		
4		

5	Intraoral Scanner Model: Infinity 5 QTY: 1 UNIT G.W(KG): Dimension(CM): X X  TPC Advanced Technology  851 S. Lawson St. Los Angeles, CA 91748  MedNet EC-REP GmbH Borkstrasse 10, 48163 Münster, Germany   	
6		
7	PI966-II	

3 Internal packaging



TPC



4 Packing box






5 Transportation Packaging




6 Shape and Structure - Dimensions

6.1 Scanner Handpiece (Main Unit)




TPC



Size	Lengths: 269mm (When installing the Scanner Tip) 204mm (Without Scanner Tip installed) Width: 39mm Thickness: 34mm
Weight	249g


6.2 Base




Size	Lengths: 129mm (Including the longest part) Width: 93mm (Including the widest part) Thickness: 149mm
Weight	258g

6.3 Scanner Tip

A) Standard




B) Mini




Size	Standard: Lengths: 116mm (longest part) Width: 33mm (widest part), Thickness: 39mm (Tail) , 17mm Mini: Lengths: 115mm (longest part), Width: 33mm (widest part) Thickness: 39 (Tail) 12.8mm
Weight	Standard: 21.4g Mini: 18.8g

6.4 Protective Tip



Size	Lengths: 87mm Width: 33mm (widest part), Thickness: 39mm (Tail) ,
Weight	16.7g

6.5 Charing Holder



Size	Lengths: 140mm Width: 114mm Thickness: 40mm
Weight	187g

6.6 Battery



Size	Lengths: 82.3mm Width: 26mm Thickness: 24.2mm
Weight	65.2g
Specification	Nominal Voltage 标称电压 3.635V Charge Voltage 充电电压 4.2V Rated Capacity 额定容量 3400mAh

6.7 Power adapter



Size	Lengths: 87mm Width: 58mm Thickness: 30mm
Weight	158g
Specification	Input: 110-240V~50/60HZ,500mA Output:12.0V $\overline{\square}$ 1.5A

6.8 U Disk with Software



Size	Lengths: 35mm Width: 16mm Thickness: 6mm
Weight	2.3g

6.9 Wi-Fi adapter



Size	91mm× 32mm× 14mm
Weight	28.2g



6.10 Practice model



Size	58mm×50mm×38mm
Weight	29.5g


6.11 Certificate/Warranty Card



Size	Lengths: 89mm Width: 50mm
Weight	2g



6.12 Packing box




Size	Dimensions inside the box: 388 mm (Length) Dimensions inside the box: 331 mm (Width) Dimensions inside the box: 142mm (Height)
Weight	1.4kg




6.13 Transportation box



A) Transport box



B) Foam





Size	Transport box: Lengths: 418mm Width: 361mm Thickness: 176mm Foam A Lengths: 388mm Width: 331 mm Thickness: 10mm Foam B Lengths:331mm Width: 142mm Thickness: 10mm Foam C Lengths: 388mm Width: 142 mm Thickness: 10mm
Weight	936g