

TPC

Laguna 2.0

Electromechanical Patient Chair

**TPC
5 Year
Limited
Warranty**



Standard Features:

- Thin back with lumbar support
- Anatomically correct back pivot point
- Adjustable double articulating headrest with contoured cushion
- Four manual chair movements and 2 programmable positions are standard
- Arms rotate for easy left or right patient entry
- Synchronized back and toe board recline
- Upholstery is designed for easy and quick clean-up
- Heavy vinyl scuff guard for toe board
- Sturdy all-metal construction
- All pivot points have permanently lubricated bearings
- Electromechanical lift system
- Quiet and dependable operation
- Chair functions are digitally controlled by a microprocessor
- Automatic return to exit
- Instant override on automatic functions
- Chair is designed for maximum serviceability
- Safety shut-off on base
- Precision photo-electric positioning