

# TPC

# Laguna 2.0

Chair Mounted Operatory System



Model #LP2115-L550LED with touchpad



Model #LP2100 with touchpad assistance side

# Laguna 2.0 Operator

## Model L2000 2.0 Dental Chair

### Standard Features

- Two operator-programmable position, one preset exit position
- Electronic foot switch preposition / auto return switch
- Fuse / surge protection
- Toe guard protection
- Synchronized seat and toe board tilt
- All metal construction
- Wide back with lumbar support
- Double articulating headrest
- Safety shut-off plate
- Pivoting armrest
- Seamless asepsis upholstery
- Input voltage (110 V or 220 V)- specify
- 5 year warranty

## Mirage 2.0 Chair Mounted Delivery System

### Standard Features

- Asepsis automatic control for three handpieces with syringe
- Individual handpiece pressure adjustments
- Non-retracting water coolant with on/off
- Individual water coolant flow controls
- Air coolant flow control
- Handpiece tubing flush system
- Ultra soft handpiece tubing
- 4 hole straight gray asepsis tubing with Midwest connector
- Stainless steel tray with non-skid pad
- Disc foot control wet/dry toggle with 7' of smooth tubing
- Equipment utility center with asepsis housing
- Vacuum canister with solids collector screen
- Straight light post for unit mount dental light
- Pre-wiring for fiber optic light source
- Air-actuated counterbalanced arm lock

## Mirage 2.0 Post Mount Utility Center

### Standard Features

- Gravity drain swing away cuspidor
- Adjustable timed bowl rinse with push-button cup filler
- Self contained water system with city/ bottle selector
- Water outlet quick disconnect
- Vacuum canister with solids collector screen

## Assistant's Instrumentation

### Standard Features

- Rotating assistant's instrument holder
- One autoclavable HVE with asepsis tubing
- One autoclavable saliva ejector with tubing

## Dental Operatory Light

**L550-LED** LED Operatory Light, 8,000 - 30,000 lux

**L600-LED** LED Operatory Light, 3,000 - 35,000 lux

**L690-LED** LED Operatory Light (Ceiling), 3,000 - 35,000 lux

## Factory Options

	Item Number
Additional 3 way syringe	#TPC3WAY+T
Additional HVE with tubing	#P2335+TB
Power optic system, ISO B (5-hole)	#H7000
Power optic system, ISO C (6-pin)	#H7660
Built-in piezo scaler with transformer	#800N
Super thin x-ray viewer- unit mounted	#XV100/150
Standard LCD monitor bracket	#WIC800
Deluxe LCD monitor bracket	#WIC825
Mirage LCD monitor bracket	#WIC840
17" standard multi media LCD monitor	#M017
Plastic utility shelf- (grey, white or beige)-post mounted	#US100
Touch pad-unit head mounted	#TP-2005

## Upholstery Selections:



**L2000 2.0** (Dental Chair)  
**2000 2.0** (Unit w/ Cuspidor)  
**2015 2.0** (Unit w/o Cuspidor)  
**L550-LED** (LED Operatory Light)  
**L600-LED** (Sensor Operatory Light)  
**L690-LED** (No Sensor Light)

