



### Features:

- Tri axis light head rotation
- Handles are conveniently located for operator positioning and asepsis protection
- 4 LED light source for illumination
- 2 LED light source
- Built in dimmer function
- PoliBlock filter mode
- Motion sensor / On - Off toggle standard
- 2 year warranty
- Packing dimensions: 37" x 17" x 15" / 28 lbs.



### SPECIFICATIONS

Input Voltage:	110/220 V, 60/50 Hz
Output Voltage:	24 Vac
Bulb Type:	LED
Light Output:	8,000 Lux High 30,000 Lux Low
Color Temperature:	5,000° Kelvin
Focal Distance:	700mm
Light Pattern:	100mm X 180"