

TPC

Mirage

Mirage 2.0 Mobile Doctors Cart

#MC-501



#MC-500
(without Vacuum)





MC-500 Mirage 2.0 Mobile Cart
MC-501 Mirage 2.0 Mobile Cart with Assistant's Pack-

Standard Features

	MC-501	MC-500
Automatic control for three handpieces	•	•
Master ON/OFF	•	•
Individual handpiece pressure controls	•	•
Individual air and water coolant controls	•	•
Three-way syringe	•	•
Easy-clean work surface, 16 1/2" x 9 1/2"	•	•
Adjustable height mobile frame, 27" - 39"	•	•
Handpiece tubings	•	•
Disc type wet / dry foot control	•	•
Central vacuum assistant's instrumentation Central vacuum solids separator Central vacuum valve with easy-clean tubing Saliva ejector valve with tubing, HVE valve with tubing	•	
Mirage utility center, 15" x 15" x 6" Manual shutoff valves with complete junction box	•	•
Air water filters, regulators	•	•
Umbilical, 7'	•	•
Non-retraction water coolants	•	•
Handpiece flush system	•	•
9" x 13" stainless steel tray with silicone-skid pad	•	•

Packing dimensions: 21" x 22" x 30" / 51 lbs

Optional Accessories

		MC-501	MC-500
#800N	Advance 800N piezo scaler-built-in	•	•
#LED55BN	LED 55BN curing light-built-in	•	•
#P2355+TB	Additional Hve valve and tubing	•	
#H7000	Fiber optic light source complete-5 hole	•	•
#H7040	Additional fiber optic tubing-5 hole	•	•
#H7660	Fiber optic light source complete-ISO 6 pin	•	•
#H7045	Additional fiber optic tubing-ISO 6 pin	•	•
#LQS666K	High speed handpiece-push button	•	•
#ELK700	Low speed handpiece complete	•	•

