



OMC 2375CV-SL

Orthodontic Mobile Delivery Cabinet




Optional OMC 2200
Monitor Bracket



Rear compartment
4" W X 12"D X 20" H



Writing surface with
optional solid surface



Artic

OMC-ST81



White

OMC-ST82



Cream

OMC-ST85



Bamboo

OMC-ST86

Optional hard surface stone top

This self-contained orthodontic mobile delivery cabinet offers a pullout breadboard, two instrument drawers, one medium and two large storage drawers for materials and supplies that fully extend and quickly remove for easy cleaning. In addition, this cabinet provides a large open space for accessories with convenient access from the front. The oversized dual-caster wheels, makes positioning of the cabinet simple and easy on all floor surfaces, including carpet.

Package include the following:

- 3 HP Automatic selection
- 3-way Syringe
- Vacuum Package with Saliva Ejector and HVE
- Self contained water bottle system
- Variable speed disc-type foot control
- Water coolant on/off switch with flow control
- 1 Year warranty
- Self contained dual compressor

Standard Cabinet Features:

- Five easy to remove and clean full extension drawers
- One pull-out writing surface.
- High quality hardware throughout
- Oversized caster wheels with LOCKS
- Four stationary handles for easy cart positioning
- A removable counter - top plastic tray
- Self contained dual compressor

Packing Dimensions (L x W x H)	Weight
26" x 28" x 37"	200 lbs.

