



Rapid 5 Intraoral Scanner Package

Contents Guide





Contents

1 Components	3
2 Label	4
3 Internal packaging.....	6
4 Packing box	6
5 Transportation Packaging	7
6 Shape and Structure - Dimensions.....	8
6.1 Scanner Handpiece (Main Unit)	8
6.2 Base	9
6.3 Scanner Tip	10
6.4 Protective Tip	11
6.5 Data Cable(Type-C).....	12
6.6 USB to Type-C adapter	13
6.7 U Disk with Software	14
6.8 Certificate/Warranty Card	15
6.9 Packing box	16
6.10 Transportation box	17



Product Name: Intraoral Scanner

Model: Rapid 5


Packaging unit: 1 set






1 Components

No.	Component name	Specifications	Quantity	Pictures
1	Scanner Handpiece	Rapid 5	1	
2	Scanner Tip	/	4	
3	Protective Tip	/	1	
4	Cradle	/	1	
5	Data Cable (Type-c)	/	1	
7	Certificate Warranty Card	/	1	

8	U Disk with Software	/	1	
9	USB adapter	/	1	
10	Packing box	/	1	


2 Label

No.	Label	Location
1	 Intraoral Scanner Model: Rapid 5 SN MD CE FC	

2		
3		
4		
5		



3 Internal packaging



4 Packing box






5 Transportation Packaging




6 Shape and Structure - Dimensions

6.1 Scanner Handpiece (Main Unit)



Size	Lengths: 246.0mm (When installing the Scanner Tip) 178.0mm (Without Scanner Tip installed) Width: 40mm Thickness: 34mm
weight	147g(Without Scanner Tip installed)


6.2 Base




Size	Lengths: 131mm/181mm (Including the longest part) Width: 43mm/70mm (Including the widest part) Thickness: 54mm
weight	216g

6.3 Scanner Tip

A) Standard




B) Mini



Size	Standard: Lengths: 116mm (longest part) Width: 33mm (widest part), Thickness: 39mm (Tail) , 17mm Mini: Lengths: 115mm (longest part), Width: 33mm (widest part) Thickness: 39 (Tail) 13mm
weight	Standard: 21.4g Mini: 19.6g

6.4 Protective Tip



Size	Lengths: 87mm Width: 39mm (widest part), 25mm Thickness: 33mm (thickest part), 20mm
weight	16.7g



6.5 Data Cable (Type-C)



Size	2m
weight	67g (with cable structural component)



Size	Lengths: 40mm Width: 14mm Thickness: 4mm
weight	6.60g

6.7 USB Disk with Software



Size	Lengths: 35mm Width: 16mm Thickness: 6mm
weight	2.3g




6.8 Certificate/Warranty Card



Size	Lengths: 89mm Width: 50mm
weight	2g

6.9 Packing box




Size	Dimensions inside the box: 388 mm (Length) Dimensions inside the box: 331 mm (Width) Dimensions inside the box: 142mm (Height)
Weight	1.4kg


6.10 Transportation box



A) Transport box



B) Foam



Size	Transport box: Lengths: 418mm Width: 361mm Thickness: 176mm Foam A Lengths: 388mm Width: 331 mm Thickness: 10mm Foam B Lengths: 331mm Width: 142mm Thickness: 10mm Foam C Lengths: 388mm Width: 142 mm Thickness: 10mm
Weight	936g