

— MIRAGE 2.0 HYGIENE PACKAGE · MHP1001

Mirage 2.0

Hygiene Operatory System, Steel Base

Hygiene operatory with the 4000-2.0 hydraulic steel base chair and H1001 rear radius delivery unit, with unit head height adjustment from 27" to 39".



Specifications

01 MODELS

Model	Description
MHP1001	Mirage 2.0 Hygiene Operatory System (without light)

02 PACKAGE INCLUDES

Part	Description
4000-2.0	Hydraulic patient chair with steel base, wide back and full lumbar support, rotating armrests left or right, double articulating headrest with sliding height adjustment, 30 degree swivel from left and right of center, 4 programmable positions with automatic return and foot control
H1001	Rear radius asepsis 3 Handpiece automatic control with coolant controls, diaphragm block design, wet/dry variable speed foot control, deluxe utility center with air /water filters and regulators, master shut off valves, Unit Head Height Adjustment from 27" to 39" Includes 3 way syringe and clean water system, with touch pad control

03 OPTIONS AND PARTS

Part	Description
LP-100	Mirage swing mount light post (ordered w/unit)
LP-105	Mirage swing mount monitor post with monitor mounting bracket (ordered w/unit)
L570-LED	Lustrous LED Ceiling Mounted Operatory Light
L550-LED	Lustrous LED Operatory Light
LP48	Lustrous LED Operatory Light
WIC800	Standard flat screen monitor bracket
WIC840	Mirage flat screen monitor bracket with 16" extension
M017	17" LCD multi media flat screen monitor
A800N-1	Advance 800 piezo scaler – built in
H7000	Illumination System – ISO B (5 hole) – tubing, control box, bulb, power supply
H7660	Illumination System – ISO C (6 pin) – tubing, control box, power supply
ALED55BN-A	LED 5W curing light – built in to dental unit TPCRDSBRKT TPC Radius Bracket
TPE-2011	Touch pad Control - Chair & Electric Motor System Touch pad
EMB-901	Built in Electric Motor System
AS-1101	Mirage Assistants stool
DR1102	Mirage Doctor's Stool
TPC-BRKT-2.0	Mirage chair Bracket
2000-PMU	PMU with cuspidor & cup filler

Options are not factory-installed unless specified on the purchase order. Contact TPC for pricing.