



### Features:

- Tri axis light head rotation
- The contoured handles are conveniently located for operator positioning and asepsis protection
- Halogen bulb is easily removed for replacement
- The dichroic-coated reflector filters out the infrared rays and dissipates heat allowing only cool, white light to illuminate oral cavity
- The two-way intensity switch allows for correct oral cavity illumination and a reduction in eye fatigue
- Removable plastic splash shield
- Optional: motion sensor
- 1 year warranty
- Packing dimensions: 37" x 17" x 15" / 28 lbs.

### SPECIFICATIONS

Input Voltage:	110/220 V, 60/50 Hz
Output Voltage:	12 V High, 10.5 V Low
Power Output:	55 W High, 50 W Low
Bulb Type:	Tungsten Halogen Type 12 V - 55 W
Light Output:	29,000 Lux High 17,000 Lux Low
Color Temperature:	4500° Kelvin
Focal Distance:	28"
Light Pattern:	3.6" x 7.1" @ 28"

