

Model #L-7601-LED



Features:

- 2 D light head rotation



- Handles are removable and sterilizable




- Packing dimensions
37" X 16" X 15" 18 lbs

Includes:

- 24V DC power Adaptor
- LED light head
- Left and Right Light handles
- Manual On / Off Toggle
- Adjustable light intensity
- PoliBlock Filter
- Post Mounted Only
- 2 Year Warranty
- Nylon Bushing
OD 1 7/8" ID 1 3/8"

PoliBlock:

Prevents polymerization
of resins.



Specifications:

Absorbed Power:	Max 8vA
Input voltage:	110V/220V, 60/50 Hz
Output voltage:	24 Volts AC
Light intensity (Lux):	3,000 - 35,000 Lux @ 700mm
Cooling System:	Static Heat sink (no fan)
Illumination area:	100mm x 175mm @ 700mm

