



Model# **PM25**

### Features:

- 25Khz
- Precision automatic tuning circuitry
- Linear power control for exact power setting
- Accepts: Cavitron, P, SLI, TFI, and FSI inserts
- Compact and lightweight
- Auto-ranging 110-240 VAC input
- Includes one ultrasonic scaler insert and foot control
- 1 Year Warranty

*Power You Want,  
Where You Want It!*

### Specifications:

Weight .....	6lbs.
Dimensions .....	8"(W) x 8"(D) x 3"(H)
Output Power to Tip .....	Operatory adjust approx. 10-40 watts
HandPiece Cable Length.....	7.5
Power Supply to Unit .....	110-240 VAC input, 50/60Hz
Water Flow & Supply Pressure ...	18cc - 30cc per minute flow from insert tip
Supply Pressure to Scaler .....	25-60 PSIG, city water or bottle system

