




### 3550 Rear Vac Pack Parts List



No.	Description	Part Number
1	Rear Mounting Hub	3550-BRKT-TPC
2	L. Bracket	3550-BRKT-TPC
3	R. Bracket	3550-BRKT-TPC
4	Extension Arm	3550-4
5	Mounting Plate Can	3550-05
6	Extension Arm	3550-RP-6
7	Lower Arm (Tube)	3550-LWARM
8	Upper Arm (Tube)	3550-UPARM
9	5/8" to 1/2" adaptor	3550-DA-1.5
10	Female Q.D.	53047
11	Syringe Block	52012
12	Screw	3550-12
13	Washer	3550-13
14	Stop Screw	3550-14
15	Washer	3550-15
16	Bushing	3550-16
17	Dome Washer	3550-17
18	90 Degree Barb	3550-18
19	Nylon Bushing	3550-19
20	Suction Canister	53008-R
21	Top Cover	3550-21
22	4 position holder	53041
23		
24	Allen Screw	3550-24
25	Screw	3550-25
26	Allen Screw	3550-26
27	Screw	3550-27
28	Bolt	3550-28
29	Nut	3550-29
30	Nut	3550-30
31	Spring	3550-31

TPC 851 S. Lawson St. City of Industry, CA 91748  
 Phone 626-810-4337 Fax 626-810-4245