



PC 2530

Mobile Vacuum System

Mobile Evacuation Unit

The PC 2530 Mobile Vac System is perfect where central installation is impractical.

Moves with Ease

Reinforced large polyethylene carpet casters

Efficient

Variable controlled motor speed can evacuate at the rate of 75 cfm, or may be adjusted to the user's preference

Simple Maintenance

Non-corrosive collection canister cleans easily and is equipped with an automatic overflow shut-off

Perfect Backup

Invaluable backup when the central system breaks down. The user just rolls it out, plugs it in, and continues operating.

Cost Effective

Low cost leader of oral evacuation units. Excellent for limited budgets.

Compact

Measures only 14" x 12" x 31" . Easily placed at chair side and stores inconspicuously.



Specifications

Electrical 7/8 HP
115 VAC, 50/60 Hz
230 VAC, 50/60 Hz
4.0 Amps

The PC 2530 Mobile Vac System is a portable, self contained vacuum system suitable for use by the dentist in a temporary location or when it is not feasible to install a central vacuum system.

Package include the following:

- Equipped with a muffler system for quiet operation
- Hve and SE valves on each suction hose
- Activates only when Hve or SE valves are lifted from hanger
- High volume suction
- White high impact resistant polyethylene case
- 1 year warranty



Packing Dimensions (L x W x H)	Weight
17" x 16" x 33"	70 lbs.