



- 4 GRAVITY DRAIN HAS A 1-1/2 NOMINAL PIPE PROTRUDING FROM FINISHED FLOOR
NOTE: PLACE TRAP IN LINE TO CONFORM WITH LOCAL CODES, SUPPLIED BY CONTRACTOR
 - 3 PLUMBING CENTRAL VACUUM TO THE UTILITY CENTER SHOULD BE SPECIFIED BY THE CENTRAL VACUUM SUPPLIER AND TERMINATED IN UTILITY CENTER WITH 3/4 O.D. TUBE PERPENDICULAR TO THE FLOOR, SIMILAR TO THE DRAIN CONNECTION
 - 2 MAKE SURE ALL PLUMBING CONFORMS TO PREVAILING LOCAL CODES
 - 1 LOCAL REGULATIONS PROVIDE THAT LICENSED PLUMBERS & ELECTRICIANS SHALL INSTALL UTILITIES
- NOTES: