
HEIGHT CHART

Gas cylinder height chart

Match the gas cylinder to the height of the person using the stool. A longer stroke gives the seat more travel and a higher top position.



SMALL · 80 MM STROKE

Users up to 5′4″

DR1102-3-SHRT

Optional on all stools.

MEDIUM · 140 MM STROKE

Users 5′4″ to 6′1″

DR1102-3

Standard on all doctor's stools.

TALL · 200 MM STROKE

Users 6′2″ and taller

AS1101-3

Standard on all assistant's stools. Optional on doctor's stools.