
BELLAGIO · REPLACEMENT PARTS

Bellagio delivery system parts list

Chair-mounted delivery system: doctor's unit head and flex arm, assistant arm, and assistant side box.

01 · Unit head and flex arm



How to order

Order by TPC part number. Call (800) 560-8222 or email service@tpcdental.com with the part number, model and serial number. Parts used in more than one place have one part number.

In this list

Unit head and flex arm · 35 items	Page 1–2
Assistant arm · 10 items	Page 3
Assistant side box · 36 items	Page 4–5

01 · UNIT HEAD AND FLEX ARM


Unit head and flex arm parts

NO.	DESCRIPTION	QTY	PART NO.
1	Balanced flex arm (spring arm) assembly, 710 × 75 × 58 mm	1	80101
2	Main arm trim cover, aluminum, Ø60 × 8 mm	1	80102
3	Main arm retaining block, brass, Ø60 × 38 mm	1	80103
4	Adjusting shim, brass, 20 × 6.8 × 3 mm	2	5015-RP-H4
5	Main arm, 760 × 380 mm	1	5015-RP-H5
6	Unit head arm trim cover, aluminum, Ø60 × 19 mm	1	80106
7	Unit head arm retaining block, brass, Ø60 × 48 mm	1	80108
8	Arm washer, white POM, Ø39 × Ø30.3 × 2 mm	1	80109
9	Retaining ring (C-clip), external, for Ø30 mm shaft	1	80110
10	Flow control needle valve, nickel-plated brass, M5	5	52007
11	On/off 3-way valve	1	52009
12	Unit head support arm, 420 × 230 mm, Ø60 mm tube	1	5015-RP-H12
13	Exhaust canister (clear bowl, cap and mounting nut)	1	80121
14	Handpiece hanger valve, square, normally open	4	5015-RP-H14
15	Hanger valve mounting block, aluminum	4	5015-RP-H15
16	Piezo LED scaler handpiece	1	A746
17	Advance LED-55B built-in curing light, 24 V, silver	1	ALED55BN
18	3-way syringe, stainless steel handle	1	5015-RP-H18
19	Instrument holder housing (upper shell)	1	5015-RP-H19
20	Keypad PCB with membrane overlay	1	5015-RP-H20
21	Unit head label panel, clear, 335 × 83 mm	1	5015-RP-H21
22	Pressure gauge with bracket	1	5015-RP-H22
23	Silicone tray mat, gray, 383 × 309 mm	1	5015-RP-H23
24	Brake release push button	2	5015-RP-H24
25	Brake button mounting sleeve, stainless steel, Ø25 × 18 mm	2	5015-RP-H25
26	Unit head top cover	1	5015-RP-H26
27	Unit head mounting plate, steel, 365 × 225 × 4 mm	1	5015-RP-H27
28	Syringe block	1	80210
29	Momentary flush valve	1	52008
30	Timer PCB with LCD display	1	5015-RP-H30
31	Unit head bottom cover	1	5015-RP-H31
32	Non-retraction water relay valve	1	52011
33	Handpiece control block, 4-position diaphragm (quad block)	1	5015-RP-H33
34	Air manifold, 9-port, nickel-plated brass	1	5015-RP-H34
35	Terminal block, 6-position, 15 A, with cover	1	80226

02 · ASSISTANT ARM

Assistant arm parts

Assistant instrument holder with HVE, saliva ejector, 3-way syringe and touch pad.




NO.	DESCRIPTION	QTY	PART NO.
1	Saliva ejector valve, short, aluminum	1	P2340
2	HVE valve, long, autoclavable	1	P2335
3	Assistant instrument holder top cover, white, 240 × 220 × 61 mm	1	5015-RP-A3
4	Assistant instrument holder base plate, aluminum	1	5015-RP-A4
5	Arm washer, white POM, Ø40 mm	2	5015-RP-A5
6	Assistant arm, outer section, aluminum, 400 × 40 × 120 mm	1	5015-RP-A6
7	Assistant arm, inner section, aluminum, 400 × 40 × 70 mm	1	5015-RP-A7
8	Touch pad: chair, operatory light, cup fill and bowl rinse	1	TP-2005
9	3-way syringe, stainless steel handle Same part as unit head item 18	1	5015-RP-H18
10	Silicone tray mat, gray, 184 × 131 mm	1	5015-RP-A10

03 · ASSISTANT SIDE BOX

Assistant side box

Side box with water bottle system, valves, transformer and support post.



03 · ASSISTANT SIDE BOX

Assistant side box parts

NO.	DESCRIPTION	QTY	PART NO.
1	Arm post trim ring, white, Ø78 mm	1	5015-RP-P1
2	Light post trim ring, white, Ø68 mm	1	5015-RP-P2
3	Light post cap, white, Ø57 mm	1	80202
4	Cuspidor port cap, white, Ø62 mm	2	80203
5	Side box top cover, white, 600 × 200 × 75 mm	1	5015-RP-P5
6	Side box top frame, aluminum, 400 × 95 × 55 mm	1	5015-RP-P6
7	Support post, steel, Ø75 × 521 mm	1	5015-RP-P7
8	Toroidal transformer, 115 / 230 V in, 25 V / 45 W out	1	5015-RP-P8
9	Adjustment valve bracket	1	5015-RP-P9
10	Solenoid valve, 24 V AC, 50/60 Hz	2	80209
11	Syringe block Same part as unit head item 28	1	80210
12	Side box right cover, white, 309 × 200 × 127 mm	1	5015-RP-P12
13	Solids collector canister, 2⅛", side port	1	80208
14	Filter bracket	1	5015-RP-P14
15	Side box frame, aluminum, 400 × 299 × 150 mm	1	5015-RP-P15
16	Water inlet quick disconnect, female, ¼"	1	53047
17	Stop plate, steel, 80 × 20 × 3 mm	1	80211
18	Magnet plate, steel, Ø20 × 7 mm	4	5015-RP-P18
19	Magnet, Ø12 × 5 mm, M4 countersunk	8	5015-RP-P19
20	Side box front / back panel, white, 389 × 311 mm	2	5015-RP-P20
21	Clean water bottle, 1000 mL	1	80218
22	Pressure gauge, 0–150 psi / 0–1 MPa, M5	1	80216
23	Regulator bracket	1	80214
24	Selector valve bracket	1	5015-RP-P24
25	Water bottle cap assembly	1	80219
26	Routing valve (toggle)	1	53014
27	On/off 3-way valve Same part as unit head item 11	1	52009
28	Side box left cover, white, with printed labels	1	5015-RP-P28
29	Water bottle cap bracket	1	5015-RP-P29
30	Terminal block, 6-position, 15 A, with cover Same part as unit head item 35	1	80226
31	Pressure regulator, nickel-plated brass, 1 in / 2 out	1	5015-RP-P31
32	Magnet plate, steel, Ø20 × 13 mm	4	5015-RP-P32
33	Side box control PCB with membrane overlay	1	80234
34	In-line fuse holder, 5 × 20 mm	1	80232
35	Air-electric switch	1	80230
36	Terminal block, 4-position, 15 A, with cover	1	80231