
H1001 · REPLACEMENT PARTS

H1001 hygiene delivery unit parts list

Rear-mounted radius delivery unit for the Mirage (bracket 16-1) or Laguna (bracket 16-2) chair.



How to order

Order by TPC part number. Call (800) 560-8222 or email service@tpcdental.com with the part number, model and serial number. Parts shared with the Mirage 2.0 and MC500 units keep the same part number.

Assemblies

Hygiene rear radius delivery unit, complete	H1001
Hygiene PMU with cuspidor	H1001-PMUC
Laguna 2.0 hygiene operatory system	LHP1001-2.0

PARTS TABLE

H1001 parts

NO.	DESCRIPTION	QTY	PART NO.
1	Non-skid tray mat, gray silicone, 419 × 239 mm	1	80133
2	Tray top shell, white, with printed logo, 450 × 350 × 60 mm	1	80132
3	Water bottle cap bracket, steel, 72 × 30 × 42.5 mm	1	MC500-RP-3
4	Water bottle cap assembly	1	80219
5	Tray bottom shell, white, with printed logo, 450 × 350 × 55 mm	1	MC500-RP-5
6	Handle end cap, stainless steel, Ø22 × 15 mm	1	MC500-RP-6
7	Tray handle, 7/8" OD	1	80117N
	Tray handle, 5/8" OD (original style)		80117
8	Flow control needle valve	4	52007
9	Momentary flush valve	1	52008
10	On/off 3-way valve	2	52009
11	Clean water bottle, 1000 mL	1	80218
12	L-shaped support arm, steel, 590 × 330 × 85 mm	1	H1001-RP-12
13	Connecting plate, 15", steel, 382 × 85 × 16 mm	1	H1001-RP-13
14	Pivot nut, steel, Ø70 × 12 mm, M45 × 2	2	H1001-RP-14
15	Thrust needle bearing, AXK4565, 45 × 65 × 5 mm	4	H1001-RP-15
16-1	Chair mounting bracket, Mirage, steel, 272 × 190 × 44 mm	1	H1001-RP-16-1
16-2	Chair mounting bracket, Laguna, steel, 200 × 190 × 83 mm	1	H1001-RP-16-2
17	Pivot shaft, steel, Ø85 × 70 mm	2	H1001-RP-17
18	Star knob screw, M8 × 25 mm, black, Ø50 mm head	1	MC500-RP-15
19	Post sleeve, white POM, Ø48 × Ø42 × 59.5 mm	1	H1001-RP-19
20	Telescoping post, stainless steel, Ø42 × 240 mm	1	H1001-RP-20
21	Post sleeve, white POM, Ø60 × Ø42 × 70 mm	1	H1001-RP-21
22	Mounting sleeve, aluminum, Ø60 × Ø48 × 80 mm	1	H1001-RP-22
23	Handpiece holder with hanger valve, normally closed	3	MC500-RP-29
24	Syringe holder, no hanger valve	3	MC500-RP-18
25	Handpiece holder bar, through-screw type (current)	1	80122N
	Handpiece holder bar, set-screw type (earlier units)		80122
26	Holder bar bracket, aluminum, 150 × 40 × 25 mm	1	80127
27	Syringe block	1	80210
28	Pressure gauge, square, front mount	1	80128
29	Routing valve (toggle)	1	53014
30	Touch pad: chair, operatory light, cup fill and bowl rinse	1	TP-2005
31	Filter bracket, steel, 60 × 60 × 7 mm	1	MC500-RP-24
32	Solids collector canister, 2⅛", side port	1	80208
33	Non-retraction water relay valve	1	52011
34	Check valve	4	52019
35	Handpiece control block, 3-position diaphragm (tri block)	1	80136-TRI
	Tri block diaphragm		80136-DIAPH-3
36	Exhaust canister (clear bowl, cap and mounting nut)	1	80121