
L700-LED · REPLACEMENT PARTS

L700-LED light head parts list

Vivid L910 LED light head used on the L700-LED, L795-LED and L700-LED-RAD operatory lights.



How to order

Order by TPC part number. Call (800) 560-8222 or email service@tpcdental.com with the part number, model and serial number.

Assemblies

Vivid LED light head, complete	L910
Vivid LED operatory light, post mount	L700-LED
Vivid LED operatory light, ceiling mount	L795-LED
Vivid LED Radiant light, post mount	L700-LED-RAD
Light head adapter to Faro light arm	L570-LHA2FA

PARTS TABLE

L910 light head parts

NO.	DESCRIPTION	QTY	PART NO.
1	Yoke, white, 248 × 240 × 48 mm	1	L700-LED-RP1
2	Yoke silicone strip, gray	1	L700-LED-RP2
3	Touch keypad PCB	1	L700-LED-RP3
4	Back housing, white, with engraved buttons	1	L700-LED-RP4
5	Phillips screw, M4 × 10 mm	2	L700-LED-RP5
6	Heat sink	1	L700-LED-RP6
7	Phillips screw, M4 × 16 mm	4	L700-LED-RP7
8	Optical module assembly (high-brightness LEDs)	1	L700-LED-RP8
9	Front housing, white, with printed logo	1	L700-LED-RP9
10	Phillips screw, M3 × 8 mm	4	L700-LED-RP10
11	Trim panel, top and bottom	2	L700-LED-RP11
12	Front housing silicone plug	4	L700-LED-RP12
13	Trim panel, left and right	2	L700-LED-RP13
14	Plastic screw, PT4 × 10 mm	4	L700-LED-RP14
15	Infrared (proximity) sensor assembly	1	L700-LED-RP15
16	Toggle switch, M6	1	L700-LED-RP16
17	Yoke pivot assembly	2	L700-LED-RP17
18	Handle shell, gray, 144 × 39 × 20 mm	2	L700-LED-RP18
19	Handle core, gray, 133 × 39 × 29 mm	2	L700-LED-RP19
20	Light head main PCB	1	L700-LED-RP20
21	Power cord, 2-wire, 60 cm	1	L700-LED-RP21