
MC500 / MC501 · REPLACEMENT PARTS

MC500 / MC501 parts list

Mirage mobile delivery cart with unit head, telescoping post and H-base on locking casters.



How to order

Order by TPC part number. Call (800) 560-8222 or email service@tpcdental.com with the part number, model and serial number.

Shared parts

Parts this cart shares with the Mirage 2.0 swing mount and Bellagio units keep the same TPC part number on every list.

Umbilical

MC500, 10 ft, installed	**MC500-UMB-10**
MC500, each added ft	**MC500-UMB-FT**
MC501, 10 ft, installed	**MC501-UMB-10**
MC501, each added ft	**MC501-UMB-FT**

PARTS TABLE

MC500 / MC501 parts

NO.	DESCRIPTION	QTY	PART NO.
1	Silicone tray mat, gray, 419 × 239 mm	1	80133
2	Unit head top shell, white, with printed logo	1	80132
3	Water bottle cap bracket, steel, 72 × 30 × 42.5 mm	1	MC500-RP-3
4	Water bottle cap assembly	1	80219
5	Unit head bottom shell, white, with printed logo	1	MC500-RP-5
6	Handle end cap, stainless steel, Ø22 × 15 mm	1	MC500-RP-6
7	Unit head handle, 7/8" OD	1	80117N
	Unit head handle, 5/8" OD (original style)		80117
8	Flow control needle valve	4	52007
9	Momentary flush valve	1	52008
10	On/off 3-way valve	2	52009
11	Clean water bottle, 1000 mL	1	80218
12	Frame tube clip, steel, 56 × 35 × 26 mm	1	MC500-RP-12
13	H-base, steel, 500 × 480 mm	1	MC500-RP-13
14	Caster, silent, M10 stem, with brake, gray	4	MC500-CASTER
15	Star knob screw, M8 × 25 mm, black, Ø50 mm head	1	MC500-RP-15
16	Post sleeve, white POM, Ø50 × 70 mm	1	MC500-RP-16
17	Telescoping post, aluminum and stainless steel, Ø60 × 400 mm	1	MC500-RP-17
18	Syringe holder, no hanger valve	3	MC500-RP-18
19	Handpiece holder bar, through-screw type (current)	1	80122N
	Handpiece holder bar, set-screw type (earlier units)		80122
20	Holder bar bracket, aluminum, 150 × 40 × 25 mm	1	80127
21	Syringe block	1	80210
22	Routing valve (toggle)	1	53014
23	Pressure gauge, square, front mount	1	80128
24	Filter bracket, steel, 60 × 60 × 7 mm	1	MC500-RP-24
25	Solids collector canister, 2⅛", side port	1	80208
26	Non-retraction water relay valve	1	52011
27	Check valve	4	52019
28	Handpiece control block, 3-position diaphragm (tri block)	1	80136-TRI
	Tri block diaphragm		80136-DIAPH-3
29	Handpiece holder with hanger valve, normally closed	3	MC500-RP-29