
Mirage 2.0 Chair Template

TPC

TOP VIEW



OVERALL LENGTH	CHAIR LENGTH	UTILITY EXTENSION	CHAIR WIDTH	BASE WIDTH	UTILITY BOX
74.01″ 1880 mm	**66.14″** 1680 mm	**7.87″** 200 mm	**25.20″** 640 mm	**24.60″** 625 mm	**12.59″** 320 mm

Water and air pressure gauges **¾″** O.D. PVC tube for central vacuum **1½″** nominal pipe for gravity drain (optional location for central vacuum) **Quad** electrical box

Transformer **Umbilical**

Dimensions measured from the chair centerline (blue). Utility locations are shown for planning; confirm against the local plumbing and electrical codes before rough-in. Not to scale.